

Meet Nova 2.0 — Our Next-Generation Glucose Meters

The Point of Care (POC) Team is excited to introduce **Nova 2.0**, the newest generation of glucose meters coming to all areas beginning in **March**.

WHAT'S NEW?

- **Seamless transition:** Minimal operational changes - screens and workflows will look familiar.
- **Greatly improved scanner:** Faster, more accurate scanning for a smoother workflow.
- **Enhanced functionality:**
 - **Built-in Wi-Fi** for improved connectivity.
 - **Integrated battery** — no more battery replacements needed.
 - **Side-button reset** to quickly troubleshoot or reboot the meter.

WHAT'S DIFFERENT?

IMPORTANT: Different order of sample type options


Patient Screen Selection Difference

The order and layout of the patient screen selection (arterial, capillary, venous) differ from the previous devices.

Legacy: Capillary, Venous, Arterial

NEW: Arterial, Capillary, Venous

- **Visual comparison:** See the diagram below for a side-by-side image of the new meter and the old meter.



- This screen configuration is system-defined and cannot be modified or customized.

FOR QUESTIONS OR CONCERNS:

The POC Team is here to support you during the transition.

☎ 614-685-6610

✉ path.glucose@osumc.edu

Stay tuned for more information as we approach go-live!

RESOURCES:

SharePoint Site: [Point of Care Testing Page](#)

